CO-170P Portable Single Gas Detector

NINGBO DSW ELECTRONICS CO.,LTD Room 1621 Linqiao Piazza, Ningbo City,315010, China

Phone: +86(574) 27861829 27861849 27861830 Fax: +86(574)87264906 E-mail: <u>sales@dswbrand.com</u> Website: <u>www.dswbrand.com</u>

Ningbo DSW Electronics Co., Ltd

Ningbo DSW Electronics Co., Ltd

Precautions:

To avoid personal safety injury, Instrument damage and potential dangerous accident; do not use the portable gas detector before reading this manual.

1. Brief Introduction

CO-170P portable single gas detector can make continuous detection to combustible and toxic gases. It is suitable for combustible and toxic gas leakage detection in underground pipe or mines, and keeps the workers safe, prevents the facilities from being destroyed.

The detector, adopting excellent-quality sensor, makes detection in the way of natural diffusion. It has good sensitivity and reproducibility. The detector adopts embedded MCU controller, easy to operate.

The shell adopts special high strength material and anti-smooth rubber, with the characters of water-proof and dust-proof.

2. Features and specifications

2.1 Features

Advanced MCU control with low power consumption; Low and high alarming level; Adjustable calibrating level; Excessive high concentration protection function; High concentration protection for combustible gas; Self test for the combustible gas sensor; Low battery indication; Replaceable sensor module; Self-adjustment function Visual and audible alarm with vibration; STEL and TWA alarming for toxic gases Advanced self-examination and self-renovation function



10.6 To assure the testing accuracy, the detector should be calibrated termly, and the calibration period should not more than one year.

- 10.7 Please send the abandoned Li batteries from the detectors to the appointed places or our company. Don't throw them into the dustbin at random.
- 10.8 Any malfunction not being included in this manual, please contact us for solutions.

Affix. Table1

Model	Range	L-alarm range	H-alarm range	L-alarm	H-alarm	TWA	STEL
□CO-170P	0-100%LEL	10% LEL \sim	25% LEL \sim	20%LEL	50%LEL		
-CH ₄		25%LEL	80%LEL				
□CO-170P	0-100%LEL	$10\% LEL \sim$	25%LEL \sim	20%LEL	50%LEL		
-C ₃ H ₈		25%LEL	80%LEL				
□CO-170P	0-100%LEL	10% LEL \sim	25%LEL \sim	20%LEL	50%LEL		
-H ₂		25%LEL	80%LEL				
□CO-170P	0-100ppm	5ppm \sim	15ppm \sim	10ppm	15ppm	10ppm	15ppm
$-H_2S$		15ppm	30ppm				
□CO-170P	0-1000ppm	25ppm \sim	100ppm \sim	35ppm	200ppm	35ppm	200ppm
-CO		100ppm	500ppm				
□CO-170P	0-2000ppm	25ppm \sim	100ppm \sim	35ppm	200ppm	35ppm	200ppm
-CO		100ppm	500ppm				
□CO-170P	0-30% vol	16%vol \sim	$22.5\% \mathrm{vol} \sim$	19.5%vol	23.5%vol		
-O ₂		19.5%vol	24%vol				
□CO-170P	0-100%LEL	$10\% LEL \sim$	25%LEL \sim	20%LEL	50%LEL		
-C ₂ H ₅ OH		25%LEL	80%LEL				
□CO-170P	0-100ppm	20ppm \sim	30ppm \sim	25ppm	50ppm	25ppm	35ppm
-NH ₃		30ppm	70ppm				
□CO-170P	0-20ppm	3ppm \sim	5ppm \sim	5ppm	10ppm	0.5ppm	1.0ppm
-CL ₂		10ppm	15ppm				
□CO-170P	0-100ppm	1ppm \sim	3ppm \sim	2ppm	5ppm	2ppm	5ppm
-SO ₂		3ppm	10ppm				

Declaration

To keep continued product improvement, **Hanwei** reserves the right to change design features without prior notice

www.dswbrand.com

8. Troubleshooting tips:

Malfunction	Possible Reasons	Solutions
	low voltage	Charge in time
Unable to Turn on the detector	Detector down	Contact distributor
	Circuit malfunction	Contact distributor
No reaction to the	Prolong time unfinished	Wait until prolong time finished
tested gas	Circuit malfunction	Contact distributor
Inaccurate indication	Overdue sensor	Replace the sensor
	Not calibrated long time	Calibrate in time
Wrong time display	Exhausted battery	Charge and reset the time
Wrong time display	Electromagnetic disturb	Reset the time
Zero calibration function unavailable	Too much drift	Calibrate or replace the sensor
Zero point excursion function unavailable	Too much excursion	Calibrate or replace the sensor

9. Accessories

one detector, one charger, one screwdriver, one set of portable accessories, two pieces of screws, one manual, and one piece of application disk, one quality certificate and one warranty card inside this package.

10. Notes

- 10.1 Prevent the detector from falling down high places or serious vibration.
- 10.2 When there is interferential high-concentration gas, the detector may not work normally.
- 10.3 Please operate and handle in strict accordance with the introduction, otherwise the result may be incorrect or you may destroy the detector.
- 10.4 The detector should not be stored or used under the circumstance with caustic gas (such as Cl2), or be use or stored under the other rigorous circumstances (including excessive high and

low temperature, higher humidity, electromagnetic field and strong sunlight).

10.5 After long-term use, if there is dust on the surface of the detector, please clean it lightly with clean soft cloth, instead of caustic impregnant or hard things. Otherwise, the surface of the machine may be destroyed.

Password management to avoid wrong operation; Explosive proof housing

2.2 Specifications

Range: See attached table 1.

Detecting Gas: combustible gas (CH_4 , C_3H_8 , H_2) and toxic gas, oxygen, Other rare toxic gas like ammonia , NO,PH₃,NH₃,NO₂,HCN, SO₂ etc also available, ,Can be specified by the Customer in advance.,


Alarming Level: see attached table 1. Accuracy: $\leq \pm 5\%$ F.S. Response time: T<30s Indication: LCD indicates the time and state Indication of alarm, fault and low voltage with LED, sound, vibration Working Environment: Operating temperature -40° C \sim 70 $^{\circ}$ C (for combustible gas) Operating temperature: -20°C~50°C (for toxic gas) Humidity: <90%RH non-condensing Working voltage: DC3.6V Li battery; 1200mAh Working time: Combustible gas: < 8h continuously Toxic gas: <300h continuously Charging time: 4h~6h Life of sensor: 2 years Protection grade: IP65 Weight: about 200g (including battery) Dimension: 110mm×60mm×40mm

3. Configuration and functions.

3.1 Outlook

2 outside configuration

1	Visual alarm	5	Buzzer
2	Push buttons	6	Recharging socket
3	LCD display	7	sensor
4	buckle		



5. Battery charging

When the power supply is not enough or the detector cannot work normally due to the low voltage, please charge in time. When the detector is deactivated, please connect the AC connector plugs of the charger with the 220V AC power. The detector will be activated automatically and indicates that it is being charged. When the charging is completed, the voltage ico is full and there is no change. Then please remove the plug and the detector can work normally.

Attentions: You can't activate the detector when is being charged in the status of power off. In order to avoid fire or explosion, please do not charge the detector when it is working in the spot. Please try not to charge when the detector is power on, otherwise the charge speed will be affected.

6. Use and replacement of Sensor


As the sensor of the detector is modularized, so please pay attention to the life of the sensor. When it is overdue, please replace it. Every half year you need to calibrate the sensor in order to guarantee the accuracy of the detector.

You'd better replace the sensor through the distributors or the local repairing department. If there is no distributor or repairing department, please replace the sensor after you get permission from our company. The sensor should be replaced with the same modularized sensor supplied by our company.

7. Accessories


For the user easy o bring, we equip buckle, crocodile buckle and flying rings with the detector. If you need the above accessories, you can fix them on the back shell of the detector with the screw. (figure 1) . If you mainly use the buckle, please take out the crocodile buckle (figure 2), and then fix the buckle (figure 3) .





6 IH

3.2 Display elements



Number	Function	
1	Battery	
2	Lock	
3	Audio	
4	vibration	
5	Alarming conditions	
6	Warning or the maxium level	
7	numeric value	
8	Error	
9	History record ①	
10	Clock	
11	Zero calibration	
12	Calibration point or calibration status	
13	Data transmission ①	
14	Set-up of the parameters	
15	Indication of the gas consentration units	

Attentions: ① this function is only available only when the product you bought has the function of infrared communications.

3.3 Pushbuttons

pushbutton	Description	
I To active the detector, press and hold it for 3 seconds		
	Press it to cancel the the operation;	
٢	 Calibration: When the device is deactivated, press for more than 3s. To deactive the detector, press and hold it for 3 seconds 	

(figure 1)

(figure 2)

(figure 3)

www.dswbrand.com



At this time, please cover the gas response hole with the calibration cap, open the gas valve to adjust the flux as 120mL/min. Within 30 seconds, if the detector detects the calibration gas, the detector will start up the calibration date analysis procedure by itself. And the detector will automatically adjust the conversion arithmetics till all the standard gas complete the response. Then the detector saves the best arithmetics and complete the calibration. After that it enters into the L-alarm set up.

Attentions: In this mode, please avoid touching any buttons, otherwise the accuracy will be effected badly.

If the detector indication figure can't reach half the calibration figure within 30s, or the gas concentration is beyond the maximum drift range of the sensor, the E ico will light, which means the gas is not proper for calibration or the sensor is destroyed. Please replace the gas or the sensor.

5. Set up of L-alarm

In this mode, the L-alarm level can be adjusted, and the screen displays as the following chart, with the figure flickering. Press 🙆 or 🔄 to adjust the flickering figure according to vour need, and then press 😳 to complete the set up. After that it enters into the H-alarm set up.





Attentions: Attached Table 2 is the L-alarm level range.

6. Set up of H-alarm

In this mode, the H-alarm the screen displays as


figure flickering. Press flickering figure

1.12 press to then

that the detector will automatically turn off.

Attentions: Attached Table 1 is the H-alarm level range.

LC Y



according to your need, and

complete the set up. After

Ningbo DSW Electronics Co., Ltd

. — .	I To increment the display value, press it;
O	I To check the device status, Press in and the LCD will display TEM, time, STEL & TWA levels@, maxium level.
	I Press and if for more than 3s to set up the parameters
<i>6</i> —»	I To decrement the display value, press it
	I To make zero calibration Press in and in and hold for more than 3s ①
	I Confirmation of the parameters set-up
NG)	I To enable/disable the beep and vibration, press it.

1) Password is needed when this operation does. Attentions:

2 only the detector for toxic gas has this function.

4. Operating instruction

4.1. Activating the detector



To activate the detector, press The detector begins a self test below:

- 1. Display Elements Test: The LCD displays all of the screen elements.
- 2. Alarm Function Test: The audible alarm beeps, the visual alarm flashes, the backlight activates briefly, and the detector emits one vibration.
- 3. Opens the vibration and alarming indication so as to check these functions;
- 4. Indication of version number




5. Date and Time: The LCD displays the date and time automatically in the following

order.










H-alarm


7、Indicates STEL、TWA level



Attentions: The above pictures can only appears for toxic gases

8、Self-test pass

If the detector passes the self-test, it enters into a short-time count down to warm up. After that it enters normal operating mode. The LCD displays the ambient gas reading




Attentions: If the self-test fails, the LCD will display relative information. For details, please refer to the Time Error Alarming and Memorizer Alarming in the Alarming Status. If count down occurs when there is no fault, the device will choose a proper time of 3-30s according to the different sensors.

4.2. Deactivating the detector

To deactivate the detector, press "², "key, then it displays the following information:

At this time, the buzzer gives beep sound. After 3 seconds, when it displays the following figure on the screen, loosen the " key. The detector is turned off.



Attentions: When the detector is not detecting status, press 😳 continuously till it returns to the detection mode.

4.3 Menu options

Attentions: As all the parameters of this mode may endanger the safety of the operator, so please make the operation carefully Within 10s after indicating inputting the password, if there is no operation or the password is wrong, the detector will be deactivated. So please input the right password in time.

2 Zero calibration

In this mode, the detector will display a figure it detected as shown in the following chart:



At this time, if 😳 is pressed or no operation within 1 minute, the detecor will accept the present concentration as the zero point, and then enters into the mode of calibration point

set up. And press 10, it jumps over the calibration step to enter into the mode of alarming level set up, with reference to the 5 and 6 items.

- Warning: This operation should be carried out in the clean air. Otherwise the accuracy of the detector will be affected differently according to the different gas concentration in the air. When entering into the next mode, if it displays E, which means the air is not clean, or the sensor is destroyed. So please choose another adjustment place or replace the sensor.
- 3、Calibration point adjustment
- In this mode, the screen displays a set of flickering figures. Press () or () to modify this figure. This figure is the standard sample gas level with which the operator needs to calibrate the detector, as shown in the following chart:



At this time, if Si pressed or no operation within 1 minute, the detecor will accept the present concentration as the calibration gas concentration, and then enters into the mode of calibration point set up.

Attentions: The calibration point concentration range is shown in the Table 1.


4、Calibration point adjustment

In this mode, the detector displays a concentration value detected as shown in the following chart:

In the clean air, if the detecting value is not 0, choose this function to make the zero calibration.

In the mode of detection, press both input password', then release the two buttons. After that input the password to calibrate the zero point.

The detector will displays the following information in turns:



Warning: This operation should be carried out in the clean air. Otherwise the accuracy of the detector will be affected differently according to the different gas concentration in the air.

4.7 Calibration and alarm levels adjustment

- If you need to re-calibrate the detector or adjust the alarm levels, please operate by the following steps:
- 1. Enter into the calibration and alarming level adjustment mode

When the device is close, press in and in the following chart:



Only the password is correct can the detector enters into the zero calibration mode.

The menu options include:

- 1. Setup of date and time
- 2. On/off of the vibration
- 3. On/of of the key sound
- 4. Mode of communications ①
- 5. Setup of the password

Attentions: (1) this function can be available only when the product you bought has the function of infrared communications.

In the mode of normal detection, press and hold (2) and (2) simultaneously until the screen displays as following figure for 1 second, and then release the buttons. it enters into the menu.



Press 🙆 or 🔄 to choose the option you need. The following is the indication picture of each option.

display	Meanings
2005	Press 🔅 to adjust the time of the device
โร <u>อปไ</u> ่ไ	Press 🔅 to turn on or turn off the key sound
ริหลีย ู	Press 🔅 to turn on of turn off the vibration
[[ח]]	Press to enter the mode of infrared communications and transmit the data to computer
PRSS	Press to modify the password (the initial : 0000)

Attentions: ① this function can be available only when the product you bought has the

www.dswbrand.com

function of infrared communications.

After entering into each item, press $\textcircled{0}$ or $\textcircled{1}$ to change options, press $\textcircled{0}$ to confirm,
press 🖑 to exit without saving. The following is the meaning of each options:

Indication pictures	Meanings
آ ۲۴م [activate the function
	Deactivate the function

4.4 Alarms

The following table describes detector alarms and shows how the LCD looks for each alarm:

Alarming type	Indicating information
Low-alarm: I Slow tone I flash I Vibration ①	
High-alarm:	HS- 38:4/ CC
I flash	ייטקב ו
I Vibration ①	
High concentration protection: ②	
I Slow tone I flash	329"
Sensor Fault:	
I fast tone	Koff

075	L alarm:		
		SI (8) 4 (02)	
	Slow tone	الم لم 🔰	
	flash Vibration	<u> </u>	
-	, and g		
1 1 1	A alarm:		
1.	Slow tone flash	-{-{	
		• • •	
-	<u> </u>		
	r F.S.alarm:		
	Slow tone		
1	flash		
Time	error:		
men acco turn	f. If succeeded, after start-up, it enters into the u of time set up. Please set up the time ording to the local time. If failed, the device will off. Please contact the supplier for repairment.	E	
When detecting, the device will try to renovate it by itself. If failed, it will turn off. Please contact the supplier for repairment.			
 Low voltage Modified tone alarming sound per second If the serious low voltage, the ico will flicker. Then the device can continue to work for 15min at most. Please recharge it at a safe place, otherwise it will turn off. 		C	
Atte	Attentions: ① It can vibrate only when the vibration item is open.		
	② It is only available for combustible gas.		
Whe	When alarming continuously, you can press 😳 to turn off the audible alarm and		
vibration, while the 📢 and 👪 are flashing.			
	4.5. To check the device status, Press and , the LCD will display TEM, time, STEL & TWA levels (), maxium level		
Atte	Attentions: ① It is only available for combustible gas.		

4.6 Automatic zero calibration